

V-Cide® Chemical Vapor Sterilant Solution



Description:

V-Cide is a liquid chemical sterilant solution for use in chemical vapor sterilizers.

Reduced Odor:

- Improved formula reduces odor-causing agents.

Gentle and Non-Corrosive:

- Compatible for use in Chemiclave® sterilizers.
- Will not corrode, rust or dull instruments when used in Chemiclave® sterilizers according to manufacturer's instructions.

Economical:

- Cost-effective for your Chemiclave®.

Effective:

- Tested with FDA clearance to market as a sterilant solution for Chemiclave® sterilizers.

Compatible:

- Safe for use with materials recommended by manufacturers for sterilization in chemical vapor sterilizers.

Intended Use:

V-Cide is a sterilant solution designed for use in chemical vapor sterilizers including the Harvey® Chemiclave®. Chemical vapor sterilizers are designed and validated to ensure sterilization using only this special solution.

V-Cide is a chemical sterilant formulated for use in Harvey® Chemiclave® sterilizers.



Product Specifications

Active Ingredient	Formaldehyde, 0.23%; Ethanol, 57.02%
Solution Color	Clear
Odor	Alcohol odor
Shelf Life	2 years

Description	Case Qty	Order#
Liter Bottle	16	VC338-1

Product Capabilities:

V-Cide significantly reduces the risk of corrosion and dull edges on instruments. V-Cide is safe for use with materials recommended by manufacturers for sterilization in chemical vapor sterilizers.

Use Precautions:

Dry instruments thoroughly before packaging and sterilization in chemical vapor sterilizer. (Residual water will interfere with sterilization and may contribute to corrosion.) Use personal protective equipment (PPE) as directed on the product Safety Data Sheet: wear eye protection, chemical resistant gloves, mask and appropriate clothing to prevent repeated or prolonged skin contact. Avoid spilling V-Cide on exposed skin or splashing in eyes. V-Cide is flammable and should be stored away from sterilizers in a dry, cool area. Keep away from open flames. Use in a well-ventilated area. When running and opening the chemical vapor sterilizer at end of cycle, avoid direct inhalation of vapors. Keep container closed between uses.

Availability:

Stocked by MicroCare Medical distributors. Contact MicroCare Medical for a distributor near you or visit MicroCare.com/Certol.

MicroCare Medical
6120 East 58th Avenue
Commerce City, CO 80022
(800) 843-3343



Visit Us:
MicroCare.com/Certol