

SpecClean™

PRE-TREATMENT FOAM SPRAY

Enzymatic



Enzymatic Pre-Treatment Foam Spray

Spec Clean™ Pre-Treatment Foam Spray is a ready-to-use neutral pH instrument pre-cleaning spray with four soil specific enzymes: protease for protein soils, amylase for starch soils, lipase for fat soils, and cellulase for fiber soils. The detergent formula also includes surfactants and an advanced corrosion inhibition system. It comes in a trigger spray bottle with a unique applicator designed to create thick long-lasting foam.

Intended Use

Spec Clean Pre-Treatment Foam Spray is applied by spraying directly on soiled instruments to start the cleaning process. It prevents soils from drying on instruments and should be applied directly to soiled items as soon as possible after use. It may also be used when there is a delay in starting the normal instrument cleaning process for spot treatment. The powerful enzyme formula initiates the cleaning process and improves worker safety by reducing the need for manual scrubbing.

Why Use Spec Clean™ Pre-Treatment Foam Spray?

- Foam stays visible, ensuring instruments are actively treated during transport or extended holding times.
- Clear gels make it difficult to determine if instruments are sufficiently covered and protected.



Multiple Applications:

- Keep soils moist during high volume delays
- Spot treat tough soils from orthopedic, gastrointestinal and laparoscopic cases
- Extra help on tough to clean instruments including serrations and box locks

- Foam never dries hard or sticky.
- Sustained 72-Hour cleaning power with four soil specific enzymes for proteins, starches, fats and fiber.
- Neutral pH and anti-corrosives make it safe for all types of instruments.

Product Capabilities

Clinical and lab testing using blood surrogates, objective soil tests and standardized surface materials including carbon steel prove:

1. Spec Clean Pre-Treatment Foam Spray dissolves soils making it easier to clean hard to reach surfaces where biofilm can form.
2. Maintains moisture on instrumentation to prevent dried soils up to 72 hours.
3. Protects instruments and prevents corrosion during extended holding times up to 72 hours.

Spec Clean Pre-Treatment Foam Spray is a viscous thick solution. It does not aerosolize when applied with foaming sprayer according to label directions within 6 to 10 inches of instruments.

Use Precautions

Since Spec Clean Pre-Treatment Foam Spray contains high foaming surfactants, instruments should be rinsed and/or flushed with coldwater to remove surfactant residue prior to further cleaning in automated washers or ultrasonic cleaning systems.

Spec Clean Pre-Treatment Foam Spray is a mild irritant and is not corrosive or caustic. At a minimum wear eye protection and non-latex gloves when applying this product at point of use locations.

Availability

Stocked by MicroCare Medical distributors. Contact MicroCare Medical for a distributor near you or visit MicroCare.com

Product Specifications

| | |
|----------------|----------------------------------|
| Dilution Ratio | Ready-to-Use |
| pH | 7.5 - 8.8 |
| Odor | Clean, fresh scent |
| Shelf life | 2 years from date of manufacture |



Spec Clean Pre-Treatment Foam Spray is a Ready-to-Use Enzymatic Spray tested for 72-Hour Cleaning Power and provides visible evidence of compliance to prevent dried soils.

Package Sizes

| Size | Case Qty | Order # |
|---------------------|----------|-------------|
| 24 oz. Spray Bottle | 15 | VPREZF240-1 |
| 1 Gallon | 4 | VPREZF128-1 |

Rev. VPREF/TDS 23240

© 2023 MicroCare, LLC. All Rights Reserved. "MicroCare Medical", "Spec Clean" are trademarks or registered trademarks of MicroCare, LLC. The information set forth herein is based on data believed to be reliable. MicroCare Medical makes no warranties express or implied as to its accuracy and assumes no liability arising out of its use by others. This publication is not to be taken as a license to operate under, nor infringe upon, any patents not herein expressly described.



6120 East 58th Avenue, Commerce City, CO 80022
(800) 843-3343 • www.microcare.com