

Spec Clean™ Endoscope Sponge

For Immediate Precleaning of Endoscopes
and Other Devices

MANUFACTURER'S INSTRUCTIONS FOR USE



Use appropriate protective equipment including non-latex gloves and eyewear while conducting cleaning procedures using Spec Clean™ Endoscope Sponge. Perform precleaning process as soon as possible after procedure.

How to use dry sponge: Open package at tear notch, remove sponge and moisten with water or detergent solution. Use sponge to aid in removal of soils from smooth surfaces of devices such as distal end of endoscope tubing, transducers, ultrasound wands or similar devices. Sponge may be used to cover distal end of insertion tube during contaminated transport. **Sponge is single use.** Do not use sponge to clean multiple items to avoid transferring soils to other devices. **Sponge is not intended to clean serrated or rough surfaces.**

How to use presaturated sponge: Open package at tear notch, remove sponge and insert in container of fresh water (up to 500 mL). Squeeze 3 to 5 times to activate fresh enzymatic solution. Squeeze sponge slightly then use moist sponge to wipe smooth surfaces of devices such as distal end of endoscope tubing, transducers, ultrasound wands or similar devices. Sponge may be used to cover distal end of insertion tube during contaminated transport. **Enzymatic solution and sponge are single use.** Do not use sponge or enzymatic solution to clean multiple items to avoid transferring soils to other devices. **Sponge is not intended to clean serrated or rough surfaces.**

Note: Presaturated sponge contains ProEZ AW Quad™ enzymatic detergent. This product is uniquely validated for enzymatic action using water at room temperature and from 65°F – 100°F / 18°C – 38°C.

How to use freshly activated enzymatic solution: Place distal end of cannulated instrument or endoscope into solution and suction through ports and channels according to device instructions. Solution is single use for flushing one device and should then be discarded.

