

SpeClean™

INSTRUMENT LUBRICANT

Ready-to-Use Spray



Ready-to-Use Instrument Lubricant

ProLube RTU is a ready-to-use spray instrument lubricant formulated to improve the functionality of jointed and hinged instruments. It is water soluble, mineral oil-free and silicone-free.

Meets AAMI ST79 Standards for Lubricants

- Water soluble and mineral oil-free.
- Completely permeable to steam and ethylene oxide sterilants.
- Non-toxic, silicone-free, food grade lubricant system.

The Clear Choice

- Water soluble formula is compatible with hot rinse cycles and steam sterilization, avoiding sticky build-up in box locks.
- Free of coloring agents that may lead to trace contamination on clean instruments before sterilization.

Easy to Use and Reduces Cross-Contamination

- Ready to use spray, no diluting necessary.
- Eliminates messy “dip” pans.



Helping to ensure the safety and well-being of your beloved animal patients.

Protects Instruments:

- Formula is a true solution and does not separate into component layers. No shaking required prior to application!
- Maintains smooth function of jointed and hinged items.
- Contains anti-corrosive and rust inhibiting agents to protect valuable instruments.

Product Capabilities

Spec Clean™ Instrument Lubricant will easily penetrate hinges and jointed areas to improve function and prevent corrosion. This product is non-toxic and does not require rinsing after application. It is water soluble and completely permeable to steam and ethylene oxide sterilants. Because the formula is not affected by high heat sterilization, it will not create sticky build-up or yellow film inside box locks and hinges. It is safe for use in all sterilizers and with all surgical grade instruments.

Spec Clean™ Instrument Lubricant eliminates the messy “milk bath” and reduces overall instrument processing time. It is quickly applied by spraying on instruments after cleaning/rinsing and before sterilization. The spray lubricant may be applied over the entire instrument surface or targeted to jointed and hinged areas. No rinsing is required after application. Allow excess solution to drain off or lightly wipe off excess.

Spec Clean™ Instrument Lubricant meets AAMI ST79 (2010-2013) section 7.5.6: “Water-soluble instrument lubricants specifically designed for compatibility with sterilization maybe used; the manufacturer’s written IFU should be followed. Instrument lubricants containing mineral oil or other oil bases should not be used, except to lubricate the internal mechanism of powered instruments as specified by the manufacturer. Oil-based instrument lubricants should not be used, because the oil will coat bacteria and the instrument surface, interfering with steam contact during sterilization.”

Package Sizes

Size	Case Qty	Order #
32 oz Spray Bottle	15	VPLR/32-1

Availability

Stocked by MicroCare Medical distributors. Contact MicroCare Medical for a distributor near you or visit MicroCare.com

Product Specifications

Active Ingredients	Non-toxic food grade lubricant with rust inhibitors and anticorrosive agents. True solution does not separate during storage
Dilution Ratios	Ready-to-use, no mixing required
pH	3.2 - 6.2
Odor	Non-specific
Color	Clear
Shelf life	3 years
Material Compatibility	Compatible with all surgical grade metals including stainless steel, titanium, aluminum, brass and plastic instruments



Apply manually after cleaning and before packaging and sterilization to improve functionality of jointed and hinged instruments.

Rev. VPLR/TDS 23277

© 2023 MicroCare, LLC. All Rights Reserved. “MicroCare Medical”, “Spec Clean” are trademarks or registered trademarks of MicroCare, LLC. The information set forth herein is based on data believed to be reliable. MicroCare Medical makes no warranties express or implied as to its accuracy and assumes no liability arising out of its use by others. This publication is not to be taken as a license to operate under, nor infringe upon, any patents not herein expressly described.



6120 East 58th Avenue, Commerce City, CO 80022
(800) 843-3343 • www.microcare.com