



SpecClean™

Kickstart Your Cleaning Process

Save Time & Money with Pre-Treatment Gel Spray

- Spray Instruments Immediately After Procedures to Keep Soils from Drying with 72-Hour Wetting Action
- Starts the Cleaning Process with Natural Enzyme Action
- Helps Avoid Spills from Water/Holding Solution During Transport



Spec Clean™ Pre-Treatment Gel Spray

Ready-to-use biodegradable spray starts soil breakdown, maintains moisture and protects instruments during transport. Neutral pH, non-corrosive formula sprays evenly and rinses off easily with no sticky residue.

Size

- 200 mL. Spray Can
- 32 oz. Spray Bottle

Order#

- VPREZG/A200-1
- VPREZG/32-1

Call Today for a Free Sample

Check out the entire line of Spec Clean™ Infection Control Products from MicroCare Medical™



PRODUCTS TO KEEP YOUR PATIENTS SAFE & HEALTHY
 Pre-Treatment • Detergents • Disinfectants • Sterilization
 800.843.3343 | MicroCare.com



Visit Us:
For product info, SDS,
TDS, Videos and more!