

MAKE BEDSIDE CLEANING QUICK & EASY

Spec Clean™ Endoscope Sponge

For Immediate Pre-Cleaning of Endoscopes and Other Devices


- Ready-to-use foam sponge is pre-saturated with correct amount of multi-enzymatic detergent to save you time.
- Spec Clean™ enzymatic detergent quickly dissolves a broad range of bioburden on external surfaces and inside channels for improved efficiency.
- Perfect for immediate precleaning of endoscopes and other medical devices including ultrasound probes and transducer wands.
- Meets professional guidelines and manufacturer instructions to wipe and flush endoscopes immediately after the device is removed from the patient.



"Precleaning should take place at the point of use immediately following the procedure. ANSI/AAMI ST91:2015"



Individually vacuum-packed, pre-saturated foam sponge is designed for single-use.



- Round cylinder, polyurethane foam sponge is contoured to fit scope surfaces.



- Spec Clean™ enzymatic detergent is a biodegradable, neutral pH and non-foaming formula.
- Non-corrosive and safe to use on all scopes.

Description	Case Qty	Order#
Round Sponge Pre-Saturated with Spec Clean™ Enzymatic Detergent 0.5 oz / 15 mL	100	QE22-1

