SAVE TIME & MONEY PROVEN TO PROVIDE 2.5X MORE COVERAGE WHILE USING LESS PRODUCT

MicroCare Medical[™] Introduces ProEZ Gel[™] Pre-Treatment Spray

This point-of-use surgical instrument pre-treatment gel helps prevent the drying of soils on dirty surgical instruments during transport to Sterile Processing.





ProEZ Gel[™] Pre-Treatment Spray has **over 2.5X wider spray pattern** than other pre-treatment gels on the market allowing for **greater and more uniform coverage**, resulting in fewer applications and **over 50% less waste of product**.







The **wider spray** also reduces the precision required to **completely cover instruments and keep soils moist** without having to manually handle the instruments.





ProEZ Gel[™]

Competitor



A unique **bio-generated enzymatic action formula** is powerful enough for every case and soil, including sticky fats. It retains moisture while preventing the deposition of 93% of organic protein soil over a 72-hour testing period.



The aerosol gel is **easy to rinse off** without leaving a sticky residue behind and the fragrance-free formula allows for a more pleasant user experience.



IMPROVE INSTRUMENT PRE-TREATMENT PROCESSES AND FULFILLS AAMI AND AORN RECOMMENDATIONS

Use ProEZ Gel[™] Pre-Treatment Spray in point-of-use treatment areas including in-patient, surgical suites, emergency departments, labor and delivery rooms, outpatient surgery centers, dental clinics, veterinary offices and other areas with reusable instruments and devices.

It is also suitable for central sterile instrument re-processing areas for pre-treatment of heavy residual soils or when high case-loads lead to cleaning delays.

ProEZ Gel[™] Pre-Treatment Spray features a unique dispenser which uses inert compressed air to deliver a wide coverage of product making it easy for surgical teams or sterile processing department staff to apply.





ProEZ Gel[™] Pre-Treatment Spray is the latest addition to the ProEZ[™] family of surgical instrument pre-treatment products.

Designed to improve productivity and protect valuable instrumentation.

Description	Case Qty	Order#
200 mL. Aerosol Spray	12	PREZG/A200-1
32 oz. Spray Bottle	15	PREZG/32-1
1 Gallon Refill	4	PREZG/128-1



Green tint for visible coverage.



PRODUCTS TO KEEP YOUR PATIENTS SAFE & HEALTHY 800.843.3343 | MicroCare.com/ProEZGel Pre-Treatment • Detergents • Disinfectants • Sterilization



